

# Spectrometer for everybody



**Pleroma** is a **spectrometer** integrated with optical MEMS chip and an image sensor developed by **Spectrochip Inc.**. Light to be measured is guided into the entrance port of Pleroma (e.g. through an optical fiber) and the spectrum measured with the built-in image sensor is output from the CSI port for data acquisition. Pleroma allow accurate measurement with low noise.

Pleroma comes with evaluation software and **Python** source code that allows settings of measurement conditions, acquiring and saving data, and displaying graphs. Pleroma is designed to take full advantage of open-source feature of **Raspberry Pi Single Board Computer (SBC)** - the most widely distributed SBC - for quick and easy customization of user applications.

# Features

- Flat-field micro concave grating chip
- Highly accurate optical characteristics
- Direct connection to Raspberry Pi SBC
- Python source code available
- Compact design for easy integration

# Applications

- Spectrophotometer
- Component analysis in food, agriculture, etc.
- Process control for chemical products
- Water content measurement
- Chromogenic reader

### > Spectral

Item		Specification	Unit
Spectral dispersion principle		Flat-field <sup>*1</sup> micro concave grating chip	-
Spectral range		300 ~ 1000	nm
Spectral resolution (FWHM)	Тур.	5	nm
	Max.	7	
Wavelength reproducibility		±0.375	nm
Stray light <sup>*2</sup>		0.04	%

\*1: Flat-field focusing is achieved by aberration-corrected concave grating design.

\*2: Stray light is measured using an image sensor with 12-bit A/D conversion.

#### > Electrical

Item	Specification	Unit
A/D conversion	8	bits
Integration time (shutter)	0~1,000,000	μs
Total noise level <sup>*3,4</sup>	1.1	%
Interface <sup>*5</sup>	CSI camera connector	-
Power consumption	158	mW

\*3: Electrical noise of this image sensor with 8-bit A/D conversion is 1.1%. Using an image sensor with more bit number of A/D conversion will greatly increase the signal-to-noise ratio.

\*4: The total noise level is the Gaussian sum of the electrical noise and the optical noise (stray light).

\*5: Connected to Raspberry Pi SBC.

### > Mechanical / Optical

Item		Specification	Unit
Dimensions (W × D × H)	Module with Stand	44.00 × 47.28 × 26.25	mm
	Module	44.00 × 35.60 × 11.60	mm
Weight (including optical fiber stand)		12	g
Image sensor type		OV9281	-
Number of spectral pixels		1,280	pixel
Input optical solid angle		24	degree
Connector for optical fiber		SMA905	-


### Temperature

ltem	Specification	Unit
Operating temperature <sup>*6</sup>	+5~+35	°C
Storage temperature	-20 ~ +70	°C


\*6: No dew condensation. When there is a temperature difference between a product and the surrounding area in high humidity environment, dew condensation may occur on the product surface. Dew condensation on the product may cause deterioration in characteristics and reliability.

# Optical Component Layout


<u>Pleroma Spectrometer</u> uses a micro concave grating chip, making it possible to deliver high throughput and highly accurate optical characteristics.



### Hg-Ar Lamp Spectrum




#### **Other Brands**



# > Application Example: Transmission Measurement

Light to be measured is guided into the entrance port of <u>Pleroma Spectrometer</u> through the sample and the acquired spectrum is output via the CSI port to a Raspberry Pi SBC. There are no moving parts inside the unit, so stable measurements are always obtained.



#### Tool Kit

Software, source code, user manual, and CAD files can be downloaded at the following link:

https://github.com/SpectroChip/spectrochip.git


If you need any support, please email to service@spectrochips.com.

Spectrochip			- 0	×
Imago Format Divitor Anolog Gain Olgital Gain ROI : X0 1001 [1]	90	Sonar Bput Fusion Nodal : 00/0281 2010 (Martin : 528) Nod 2010 (Harper II (Martin ) Nod Ster I Jun x 3 um	Sare As Defails : Soli Soli SAVE SAVE	
RCE Scen Acts Scen Internity Max Theosed (%) The billion	With Res  Starting    EXAT  Starting    2007  Date Res    201  Date Res    9  21(4 - 2)(4)    1  4532 terminis	x1  0  p  0    x2  0  p  0    x3  0  p  0    x4  0  p  0    x5  0  p  0    x6  0  p  0    x7  0  p  0    x8  0  p  0		
START ontruno 6 AJTO 7080 300 0.5 0.5 0.5 0.6 0.5 0.5 0.6 0.5 0.6 0.5 0.5 0.6 0.5 0.5 0.5 0.6 0.5 0.5 0.5 0.5 0.5 0.5 0.5 0.5	61 62 53 54 63 64 64 64		4 67 54 54	
DONE	Pod	#aoderyh		j.

### Dimension (unit: mm)

Pleroma with SMA905 Connector

#### Pleroma Spectrometer



Information described in this material is current as of December 2022.

Product specifications are subject to change without prior notice due to improvements or other reasons. This document has been carefully prepared and the information contained is believed to be accurate. In rare cases, however, there may be inaccuracies such as text errors. Before using these products, always contact us for the delivery specification sheet to check the latest specifications. The product warranty is valid for one year after delivery and is limited to product repair or replacement for defects discovered and reported to us within that one year period. However, even if within the warranty period we accept absolutely no liability for any loss caused by natural disasters or improper product use. Copying or reprinting the contents described in this material in whole or in part is prohibited without our prior permission.



Website : http://www.spectrochips.com/ Address : No. 951, Fuxing Rd., Zhubei City, Hsinchu County, 302056, Taiwan (R.O.C.) Tel : +886 3 5520892 E-Mail : service@spectrochips.com Download : https://github.com/SpectroChip/spectrochip.git